


Best Seller

2~8°C Pharmacy Refrigerator

High-quality Product Temperature Uniformity
 Ensuring Safety of Storage



YC-1505L



Temp. Uniformity



Pulling Down to 5°C



Temp. Recovery After Door Opened for 1 Min

**Typical freezer data based on internal testing with freezer set point at 5°C and ambient temperature at 25°C.*

● Features



Temperature Control System
 Display temp.in 0.1°C increments



Double-pane Glass Door
 Automative electric heating function



Hydrocarbon Refrigerant
 Enviroment-friendly Natural Gas



CE Certificate






Built-in Data Logger
 Export up to 12 months of temperature data



Alarms System
 Multiple alarm functions

● Structure





2~8°C Pharmacy Refrigerator

Model	YC-1505L
Capacity	1505L/53.15Cu.Ft
Internal Size (W*D*H) mm/Inch	1685*670*1514/66.34*26.38*59.61
External Size (W*D*H) mm/Inch	1798*832*1997/70.79*32.76*78.62
Package Size (W*D*H) mm/Inch	1928*988*2165/75.91*38.90*85.24
NW(Kgs/Lbs)	322/710.00
GW(Kgs/Lbs)	414/912.71

Performance

Temperature Range	2~8°C
Ambient Temperature	16~32°C
Cooling Performance	5°C
Climate Class	N
Controller	Microprocessor
Display	Digital Display

Refrigeration

Compressor	1pc
Cooling Method	Air Cooling
Defrost Mode	Automatical
Refrigerant	R290
Insulation Thickness(mm)	D:55, L/R:55, B:55, U/D:65/62

Construction

External Material	PCM
Inner Material	Stainless Steel
Shelves	18
Access Test Port	1pc. Ø 23mm
Casters	6(3 brakes)
Data Logging/Export	Record every 10 mins/Max.12 months via USB port
Interface	Remote Alarm Contact, RS485
Backup Battery	Yes

Alarms

Temperature	High/Low Temperature, High Ambient Temperature, Condenser High Temperature Alarm
Electrical	Power Failure, Low Battery
System	Data Communication Failure, Door Ajar, Data Logger Failure, Sensor Fault Alarm

Electrical

Power Supply(V/Hz)	220~240V;50/60Hz	115V~/60Hz
Rated Current(A)	5.2A	9.5A

*Zhongke Meiling reserves the right to change products and specifications without prior notice.