



# Best Seller

## 2~8°C Pharmacy Refrigerator

Intelligent Control & Constant Temperature  
Multiple Alarms Bring More Security



Temp. Uniformity



Pulling Down to 5°C



Temp. Recovery After Door Opened for 1 Min

*\*Typical freezer data based on internal testing with freezer set point at 5°C and ambient temperature at 22°C.*

### YC-55EL

## Features

**Temperature Control System**  
Display temp. in 0.1°C increments

**Compact and Fit**  
Fits under standard 34" high ADA compliant countertop

**HC Nature Hydrocarbon Refrigerant**  
Environment-friendly Natural Gas



**CE Certificate**

**Foaming Door**

**Alarms System**

## Structure





2~8°C Pharmacy Refrigerator		
Model	YC-55EL	
Capacity	55L/1.94Cu.Ft	
Internal Size (W*D*H) mm/Inch	444*440*404/17.48*17.32*15.90	
External Size (W*D*H) mm/Inch	540*560*632/21.20*22.00*24.90	
Package Size (W*D*H) mm/Inch	575*617*682/22.64*24.30*26.85	
NW(Kgs/Lbs)	35/77.16	
GW(Kgs/Lbs)	41/90.39	
Performance		
Temperature Range	2~8°C	
Ambient Temperature	16~32°C	
Cooling Performance	5°C	
Climate Class	N	
Controller	Microprocessor	
Display	Digital Display	
Refrigeration		
Compressor	1pc	
Cooling Method	Forcedair Cooling	
Defrost Mode	Automatical	
Refrigerant	R600a	
Construction		
External Material	PCM	
Inner Material	Sprayed Aluminum Plate	
Shelves	2(Plastic Coated Steel Wire Shelves)	
Access Test Port	1pc. Ø 23 mm	
Casters	2+(2 Leveling Feet)	
Data Logging/Export	Record every 10 mins/Max.12 months via USB port	
Interface	Remote Alarm Contact, RS485	
Backup Battery	Yes	
Alarms		
Temperature	High/Low Temperature, High Ambient Temperature	
Electrical	Power Failure, Low Battery, Sensor Error	
System	Data Communication Failure, Door Ajar, Data Logger Failure	
Electrical		
Power Supply(V/Hz)	110V/60Hz	220~240/50Hz
Rated Current(A)	0.79A	0.9A

\* Zhongke Meiling reserves the right to change products and specifications without prior notice.