




High-quality product
Seven sensors
Temperature uniformity
Ensuring safety of storage



Temp. Uniformity



Pulling down to 5°C



Temp. recovery after door opened for 1 min

*Typical freezer data based on internal testing with freezer set point at 5°C and ambient temperature at 25°C.

YC-1320CL

Features



Touch screen
7" LCD touch screen



Double-pane glass door
Automatic electric heating function



Hydrocarbon refrigerant
Environment-friendly natural gas



CE certificate



Built-in data logger
Export up to 10 years of temperature data





Alarms system
Multiple alarm functions

Structure



Dimensions

Specification



2~8°C Pharmacy Refrigerator

Model	YC-1320CL
Capacity	1320L
Internal Size (W*D*H) mm/Inch	1340×670×1515/52.8*26.4*59.7
External Size (W*D*H) mm/Inch	1453*826*1998/57.2*32.5*78.7
Package Size (W*D*H) mm/Inch	1588*991*2179/62.5*39.0*85.8
NW/GW (kgs)	258/350

Performance

Temperature Range	2~8°C
Ambient Temperature	16~32°C
Cooling Performance	5°C
Climate Class	N
Controller	Microprocessor
Display	7"HD IntelligentTouch Screen

Refrigeration

Compressor	1pc
Cooling Method	Air Cooling
Defrost Mode	Automatic
Refrigerant	R290
Insulation Thickness(mm)	D:55, L/R:55, B:55, U/D:65/80

Construction

External Material	PCM
Inner Material	Stainless Steel
Shelves	10
Access Test Port	1pc. Ø 23mm
Casters	4(2 brakes)
Data Logging/ Interval /Recording Time	USB/Record every 5 minutes/100,000 data
Interface	Remote Alarm Contact, RS485
Backup Battery	Yes

Alarms

Temperature	High/Low Temperature, High Ambient Temperature, Condenser High Temperature Alarm
Electrical	Power Failure, Low Battery, Battery Power Detection
System	Sensor Error, Data Communication Failure, Door Ajar

Electrical

Power Supply(V/Hz)	110V/60Hz	220-240V;50/60Hz
Rated Current(A)	10.55A	5.59A

*Zhongke Meiling reserves the right to change products and specifications without prior notice.