





YC-1320L

High-quality product
Seven sensors
Temperature uniformity
Ensuring safety of storage



Temp. Uniformity



Temp. Stability



Temp. recovery after
door opened for 1 min

*Typical freezer data based on internal testing with freezer set point at 5°C and ambient temperature at 25°C.

Features



Temperature control system
Display temp. in 0.1°C increments



Double-pane glass door
Automatic electric heating function



Hydrocarbon refrigerant
Environment-friendly natural gas



CE certificate



Built-in data logger
Export up to 12 months of temperature data





Alarms system
Multiple alarm functions

Structure



● Dimensions



● Specification

2~8℃ Pharmacy Refrigerator		
Model	YC-1320L	
Capacity	1320L	
Internal Size (W*D*H) mm/Inch	1340×670×1515/52.8*26.4*59.7	
External Size (W*D*H) mm/Inch	1453*826*1998/57.2*32.5*78.7	
Package Size (W*D*H) mm/Inch	1588*991*2179/62.5*39.0*85.8	
NW/GW (kgs)	258/350	
Performance		
Temperature Range	2~8℃	
Ambient Temperature	16~32℃	
Cooling Performance	5℃	
Climate Class	N	
Controller	Microprocessor	
Display	Digital Display	
Refrigeration		
Compressor	1pc	
Cooling Method	Air Cooling	
Defrost Mode	Automatical	
Refrigerant	R290	
Insulation Thickness(mm)	D:55, L/R:55, B:55, U/D:65/80	
Construction		
External Material	PCM	
Inner Material	Stainless Steel	
Shelves	10	
Access Test Port	1pc. Ø 23mm	
Casters	4(2 brakes)	
Data Logging/Export	Record every 10 mins/Max.12 months via USB port	
Interface	Remote Alarm Contact, RS485	
Backup Battery	Yes	
Alarms		
Temperature	High/Low Temperature, High Ambient Temperature, Condenser High Temperature Alarm	
Electrical	Power Failure, Low Battery, Sensor Error	
System	Data Communication Failure, Door Ajar,Data Logger Failure	
Electrical		
Power Supply(V/Hz)	110V/60Hz	220-240V;50/60Hz
Rated Current(A)	10.55A	5.59A

*Zhongke Meiling reserves the right to change products and specifications without prior notice.